



KneeScōot

Owner's Manual

www.KneeScoot.com

WHAT'S INCLUDED

- A:** Adjusting handle bar
- B:** One Attachable Storage Bag
- C:** Rear frame
- D:** Front frame
- E:** Adjustable Knee Pad
- F:** Pin
- G:** Clamp Lever
- H:** Locking Knob
- I:** Screw
- J:** Locking Knob



ASSEMBLY INSTRUCTIONS

1. Connecting the frame - Connect front frame (D) with rear frame (C) – Insert the interface of rear frame (C) to the groove of the front frame (D). Then adjust the hole position and insert the screw (I) and then tighten up the clamp lever (G) in a clockwise



2. Adjusting handle bar (A) assembly – Attach the handle bars to the steering column. Ensure the brake line is not twisted around the bar. Place the handle bars at the desired height, making sure that the holes for inserting the locking knob (H) line up. Secure the handle bars with the locking knob (H). Twist the locking knob clockwise until fully tightened.



3. Insert the knee pad post into the receiver at the top of the knee walker frame. Lock the pad in place by inserting the locking knob (J) into the desired adjustable height hole making sure that the holes for inserting the locking knob line up. Twist the knob clockwise until it is fully tightened. Then insert the pin (F) and lock.



4. Securing the brake cable - Take the brake cable and secure it inside the two Velcro straps which you will find at the rear frame (C) as well as the steering column.



OPERATING INSTRUCTIONS

General Information

Position your leg firmly and securely in the middle of the knee pad, using the entire length of the knee pad to support your leg.

Whilst your knee is on the Knee Scoot, stand as straight as possible. Adjust the height of the handlebars to a comfortable position. The height of the knee pad can also be adjusted.

Using the Brakes

To lock the hand brake simply pull the brake lever on the handle bar towards you. To lock the brake system, pull the brake lever towards you and push down on the locking pin on the top lever. This pin will lock the wheels in place.

To release the brake system, squeeze the lever to disengage the pin.

Folding the Steering Column

Undo the locking mechanism at the bottom of the steering column to allow it to be folded down. If required you can fold it down even further by remove the knee pad.



OTHER INFORMATION

SPECIFICATIONS

We are constantly adding new products and as this manual serves as a general guide for all our products, we advise you to check out www.kneescoot.com for the latest and up to date specifications.

WARNING

Always be aware of your surroundings considering any obstacles or hills.

Always lock the brakes before attempting to transfer to or from the knee walker.

Always make gentle turns on your knee walker to ensure there's no risk of imbalance

Only one person should operate the knee walker at a time.

The maximum weight for our current range of products is 136kg / 21 stone.



EC CERTIFICATE OF CONFORMITY

QUESTIONS



hello@kneescoot.com



00 353 (0)1 866 3366



www.kneescoot.com